

Each cushion is given its own unique serial number - this makes the cushion uniquely identifiable and ensures traceability. The serial number is recorded in our merchandise management system and contains information about:

- ✓ Batch number
- ✓ Date of manufacture
- ✓ Cushion type /height /size/1- or multi-chamber system
- ✓ Date of dispatch
- ✓ Recipient (specialised trade / medical supply store)


Serial number is on:

- ✓ Bottom plate of the cushion (twice)
- ✓ Single carton (packing carton)
- ✓ Packing list (in case of several packages)
- ✓ Delivery note
- ✓ Invoice


Bottom Plate
CONFORM® Med. Air Cushion System



1 x label : Article designation, serial number and date of manufacture (JJJJ/MM/TT)



1x stamp: serial number and date of manufacture (JJJJ/MM/TT)



Single Box (Packing Box)



1 x label : Article designation, serial number and date of manufacture (JJJJ/MM/TT)

