

## Productlist 2022

### CONFORM<sup>®</sup> Mattress Overlay System



The system is used as an overlay on top of the standard bed mattress.

A complete mattress system consists of 4 sections (BMC 4) which are held together by connectorstrips.

Each inflatable mattress section is equipped with 2 valves to enable inflation/deflation from either side of the bed.

The use of a bed sheet on top of the mattress is mandatory. For incontinence people a waterproof cover is recommended use a waterproof cover. In case of the cremes on the patients skin, we also recommend waterproof cover since the natural rubber might soak and get porouse.




Article No.	Code	Details
11503	BMC 2	195x85 cm 2 inflatable mattress sections 2 foam sections with waterproof cover connector strips and pump
11504	BMC 3	195x85 cm 3 inflatable mattress sections 1 foam section with waterproof cover connector strips and pump
11505	BMC 4	195x85 cm 4 inflatable mattress sections connector strips and pump

**CONFORM<sup>®</sup> Medical Air Cushion Systems are Medical Products - Class 1**  
(Reg. No. DE/CA /23/051603NA001.01; No. 0037437)

#### Accessories

Article No.	Description
14756	CONFORM <sup>®</sup> XXL Box, Height, 7.5 cm
10572	CONFORM <sup>®</sup> Pump
10623	Foam section - Height 7,5 cm (85 x 48 cm - Foam made of Polyurethan)
2405	CONFORM <sup>®</sup> single mattress section - Height 7,5 cm (width 48 x depth 85 cm)
10576	Mattress Incontinence Cover
	(100x 200 cm, Fematex waterproof, vapur permeable 2 way stretch cover with nylon backing)
10579	Connector strips

**CE** CE-Declaration: The manufacturer declares, that CONFORM<sup>®</sup> Med. Air Cushion Systems are in accordance with all requirements for medical products class 1 as set out in the guideline 2017/745/EU (Appendix VIII, Chapter III, Paragraph 4.1 - Regulation 1).

 Wolff Kunststoffe GmbH  
Industriestraße 28  
69509 Mörlenbach-Germany  
VAT No.: DE 144458680  
SRN: DE-MF-0000178949

**Certified according to:**  
ISO 9001:2015 and ISO 13485:2016

