

Each cushion variant of our CONFORM® Med. air cushions (height/dimension/1- or multi-chamber) has its own UDI DI. Each cushion is given a unique serial number after passing quality control in our production department. This results in the UDI PI with the date of manufacture. This ensures that every CONFORM® Med. air cushion is traceable. The above data is recorded in our merchandise management system and provides further information on:

- ✓ Batch No.
- ✓ Date of dispatch
- ✓ Recipient (trade / medical supply store)

The labelling on the product and packaging is illustrated and described below:

Bottom plate:
CONFORM® Med. Luftkissensysteme



Label: Article, Serial No., Manufacturing Date (JJJJ/MM/TT)	Stamp: Serial No., Manufacturing Date (JJJJ/MM/TT)
Art. Nr./Art.No. SKH4545 SN #138800 2025/03/03	138800-2025-03-03
 (01)04251815702559(11)250303(21)138800	
Label: Barcode (GS1-128): UDI PI + Figures	

Single box (Packing Box)



Label: Barcode (GS1-128): UDI PI + Figures (01)04251815702559(11)250303(21)138800
Label: Article, Serial No., Manufact. Date, (JJJJ/MM/TT) UDI DI and further details like Cells, Dimension and accessories CONFORM® Med. Air Cushion SKH4545 S#138800-2025/03/03 UDI DI 04251815702559 1-chamber cushion Cells: 9x9, Height: 10 cm Dimension: width 43 x depth 43 cm incl. stretch-cover, air inflation pump and user manual